

Product

DATA SHEET

PETSPACER FOR DOGS

Medical device for administering inhaled treatments to pets with respiratory problems such as asthma or bronchitis.

FEATURES

- Adaptable silicone mask
- Unidirectional valve
- Transparent antistatic chamber
- Flow indicator
- Compatible with all types of inhalers

TECHNICAL SPECIFICATIONS

Dimensions	Length: 14,3 cm - Diameter: 5,5 cm
Volume	150 ml
Materials	
Body	(Acrylonitrile Butadiene Styrene)
Nozzle	MABS (Methacrylonitrile Acrylonitrile Butadiene Styrene Methacrylate) antistatic
Valve	TPE (Thermoplastic Elastomer)
Adapter	TPE (Thermoplastic Elastomer)
Coupling	MABS (Methacrylonitrile Acrylonitrile Butadiene Styrene Methacrylate)
Masks	
Mask sizes	Dogs (S/M) , Cats (M/L)
Latex Free	Yes



PACKAGE CONTENTS

- 1 chamber
- 1 mask (depending on model)
- Instructions for use

MODE OF USE/CLEANING



Scan the QR to see how to use it and the cleaning instructions.

Product

DATA SHEET

PETSPACER FOR CATS

Medical device for administering inhaled treatments to pets with respiratory problems such as asthma or bronchitis.

FEATURES

- Adaptable silicone mask
- Unidirectional valve
- Transparent antistatic chamber
- Flow indicator
- Compatible with all types of inhalers

TECHNICAL SPECIFICATIONS

Dimensions	Length: 14,3 cm - Diameter: 5,5 cm
Volume	150 ml
Materials	
Body	(Acrylonitrile Butadiene Styrene)
Nozzle	MABS (Methacrylonitrile Acrylonitrile Butadiene Styrene Methacrylate) antistatic
Valve	TPE (Thermoplastic Elastomer)
Adapter	TPE (Thermoplastic Elastomer)
Coupling	MABS (Methacrylonitrile Acrylonitrile Butadiene Styrene Methacrylate)
Masks	
Mask sizes	Dogs (S/M) , Cats (M/L)
Latex Free	Yes



PACKAGE CONTENTS

- 1 chamber
- 1 mask (depending on model)
- Instructions for use

MODE OF USE/CLEANING



Scan the QR to see how to use it and the cleaning instructions.